




CAUTION

1. Ensure the ear tip is securely attached before using the Smartbud.
2. To avoid injury, DO NOT insert the tip of the Smartbud further than 1 centimetre (cm) into your ear. The average adult human ear canal is roughly 2.5 cm (1 inch) deep and ends with the eardrum (refer to Figure 1, below). Seek medical advice if you have wax or debris beyond this depth.
3. DO NOT touch the eardrum itself or attempt to remove wax / debris resting on, or attached to, the eardrum.
4. DO NOT use the Smartbud to clean your ears unless you have successfully connected the device to your mobile phone through the Android / iOS App and are streaming live video footage.
5. Proceed slowly and carefully with caution. If you feel any resistance, DO NOT apply force.
6. Correct use of the Smartbud should not cause pain or discomfort.
 - If you feel any pain or discomfort, remove the device from your ear immediately.
7. This device is only to be used to remove soft wax and loose debris.
 - If you have hard wax or debris that is lodged in place, please do not attempt to remove it with the Smartbud and seek medical advice.
8. The Smartbud is only to be used by yourself - not to be controlled by someone else. Please refer to our 'How To Use' email to understand how to use Smartbud.
9. This product is not suitable for use on or by children under the age of 18. DO NOT use on pets or animals. Keep out of reach of children.

Figure 1 - Ear Anatomy



10. DO NOT use the Smartbud if you have, or have recently had, an ear infection. OR if you have, or have recently suffered, an injury to your ear or eardrum.
11. DO NOT use the Smartbud if you are prone to involuntary movements, dizziness and / or fainting.
12. DO NOT use the Smartbud whilst under the influence of alcohol or drugs. Improper use can cause injury. If problems occur, seek medical advice.
13. Unscrew the ear tip before cleaning. Wipe carefully with alcoholic cotton swabs or hot water for best results. DO NOT use the Smartbud without an ear tip, or with an unclean or damaged ear tip.
14. Please observe your surroundings before use. Stay at least 2 feet (ft) away from people, animals, or moving objects when using the Smartbud. DO NOT use the device in a moving vehicle or other unstable environments. Please DO NOT move around when using the device.
15. DO NOT place the Smartbud in a humid environment for a long time and avoid exposing to direct sunlight. DO NOT soak the Smartbud in liquids.
16. There is a built-in rechargeable lithium-ion battery. Even if you do not use the Smartbud for an extended period(s) of time, you should charge it once-per-month to avoid long-term battery damage. Please use a 5 Volt (5V) charger to recharge the device. DO NOT use the Smartbud product when it is charging.
17. DO NOT place heavy objects on top of the Smartbud.
18. It is normal for the temperature of the device to rise slightly when in use.

If you have any questions then please contact us via email at support@thesmartbud.com.




Step 1: Download the APP

iOS : Search "Smartbud" in App Store
Android: Search "Smartbud" in Google Play




Step 2: Connect to Product WiFi

Press and hold the power button to turn on the device, and then check your smartphone, connect the WiFi "Smartbud-XXXX"



Step 3: Go to Smartbud APP

Click 'GO' to begin using your Smartbud



Replace the ear tip


Please check that the ear tip is tight before using it.



The ear tip has its own thread and can be unscrewed clockwise for replacement.

For a better user experience, we recommend replacing the ear tip at least every 3 months. (This can be adjusted according to the frequency of use.)

Accessories



Indicator light



- The red light is always on when charging. When fully charged, the red light goes out.
- Blue light flashes when Product WiFi is not connected. Blue light is always on when Product WiFi is connected.
- When low battery, blue light flashes rapidly.

Product Parameters

Brand	Smartbud
Product Name	Smart visual ear cleaner
Model	bebird R1
Pixel	3.0 megapixel
Lens diameter	3.5mm
Optimum focal length	12 ~ 15mm
Depth of field	10 ~ 50mm
Image transmission rate	20fps
Working temperature	-10 ~ 45°C
Gyroscope	3axis
Battery	180mAh
Battery life	About 1 hour
Charging time	About 1 hour
Input current	5V=0.5A
Rod weight	13g
Working frequency	2.4Ghz
Network standard	IEEE 802.11b/g

Supplier details: Smartbud, Mist Productions Ltd, PO Box 6145, Westcliff-on-Sea S51 9AF, United Kingdom.

Warranty Information

When you shop with Smartbud, you have the peace of mind knowing that all our customers have a 90-Days Warranty. This warranty only covers product(s) that have been received defective or damaged. This warranty does not cover misuse, damage resulting from an accident, theft, improper storage, neglect, physical or electrical modifications, products not sold directly from us or one of our authorised retailers or resellers. Products received defective or damaged may be returned for a replacement. To qualify for the warranty claim, we require a video that would support the issue being reported. Please send via email to support@thesmartbud.com and our team will assist you further.

Disposal

The Smartbud contains a 180mAh battery and should be disposed of properly at a hazardous waste drop off facility.



Troubleshooting Guide

Problem	Solution
Start up and no response	Please try several more times according to the correct boot-up method Please ensure the device is charged
Camera does not turn on	Enter phone settings, search for 'Smartbud' APP. Enable local network, mobile data and camera permissions
APP quits quickly	Delete the existing APP and redownload
The power indicator of the Smartbud shows blue light flashing rapidly	Charge the Smartbud
WiFi connection is unstable	Restart the Smartbud. Reconnect to the Smartbud APP
The Smartbud fails to display real-time image after successfully connecting to the Smartbud APP	Disconnect the phone from WiFi, restart the product and reconnect the WiFi
If the above operation still fails to solve the problem please refer to tutorial videos located within the Smartbud APP.	
If further assistance is needed then please contact us via email at support@thesmartbud.com	