

## QuantiVirus™ Real-Time PCR Coronavirus (SARS-CoV-2) Test Clinical Evaluation

Clinical evaluation of the QuantiVirus™ Real-Time PCR Coronavirus (SARS-CoV-2) Test (CE-IVD) was conducted with contrived sputum specimens including 32 positive and 32 negative samples. 32 sputum sample were contrived with positive control template and 32 specimen were contrived only with the negative extraction control (NEC) template. The contrived clinical samples were tested blindly to generate the Positive Percentage Agreement (PPA), Negative Percentage Agreement (NPA) and overall percentage agreement (OPA) as a measurement of estimated Diagnostic Accuracy.

	Contrived Samples			
		Positive	Negative	Total
<b>QuantiVirus™ Real-Time PCR Coronavirus (SARS-CoV-2) Test</b>	Test positive	32	0	32
	Test negative	0	32	32
	Total	32	32	64

**Table: Contrived Clinical Sample Evaluation**

**The Table Shows Assay Accuracy for the QuantiVirus™ Real-Time PCR Coronavirus (SARS-CoV-2) Test**

Clinical sensitivity = 100% (95% CI, 88.6-100%)

Clinical specificity = 100% (95% CI, 88.6-100%, 95%CI)

Positive percentage agreement (PPA) =100 %

Negative percentage agreement (NPV) = 100%

Overall percentage agreement (OPV) = 100%