
# DOSING CAPSULES

For clean and easy loading.
Dosing Capsules can be prefilled,
stored in the Magazine, and inserted into the
Filling Chamber before operation.
This way, patients have their daily doses
contained and portable.

Dosing Capsules are suitable for herbs. Liquids can be used with the additional Liquid Pad.\*



Also compatible with the VOLCANO MEDIC in combination with the MEDIC VALVE Filling Chamber Reducer.

\*Liquids are only to be used with the VOLCANO MEDIC

#### **MANUFACTURER**





STORZ & BICKEL GmbH & Co. KG Rote Strasse 1 • 78532 Tuttlingen / Germany +49 7461 96970735 • info@vapormed.com www.vapormed.com



# **Medical Cannabis Inhalation**









#### INTENDED USE

The VOLCANO MEDIC and the MIGHTY MEDIC are intended for vaporizing and subsequently inhaling cannabinoids. The goal is the application of active ingredients in the human body via the alveoli for the appropriate medical indication. It is suitable for the temporary inhalative application of cannabinoids prescribed by a physician for use at home or in a hospital.

## **EFFICIENCY**

The inhalative application of cannabinoids is by far the most preferred method to ingest cannabinoids. One of the reasons is the rapid onset of effect. Both, the VOLCANO MEDIC and the MIGHTY MEDIC ensure an effective inhalation. They vaporize cannabinoids by hot air, without combustion. Therefore, no noxious combustion substances are generated. Reference values validated in studies show that approx. 1/3 of the cannabinoids inserted into the device enter the bloodstream. Since 2009, the VOLCANO MEDIC has allowed medical inhalation therapy with pure cannabinoids dissolved in alcohol and cannabinoids directly from dried hemp flowers. VAPORMED inhalers have been used in numerous medical studies and have established their reputation as certified medical inhalers for cannabinoids, licensed for Europe, Canada and other countries.\*

\*More information on www.vapormed.com

### MEDIC VALVE BALLOON

The removable Valve Balloon allows a large freedom of movement and an application independent of the basic device. Inhalation out of the Valve Balloon is comfortable due to practically no draw resistance.



20.0 cm (7.9")

# MEDICAL AND PORTABLE

The new MIGHTY MEDIC is the first medically approved portable inhaler. It uses a full hot air convection heating combined with additional conduction, ensuring an efficient vaporization from the first draw.

The mobile, battery-powered MIGHTY MEDIC vaporizer is characterized by its ease of handling. The two high-performance rechargeable batteries ensure a long period of use and short heating-up times. It also offers a Cooling Unit for pleasant inhalation:

Great performance, pocket sized.



8.0 cm (3.2")

3.0 cm (1.2")