



Since 1991

 (+202) 26850425 - 24872105

 +20 106 521 8282

 [info@Lmc-eg.com](mailto:info@Lmc-eg.com)  [www.Lmc-eg.com](http://www.Lmc-eg.com)  415 Ramsis st., Abbassia, Cairo , EGYPT

## Gymna 200 Series

Combination device for 2-channel electrotherapy, 2-pole and 4-pole.  
Diagnostics, ultrasound, laser, simultaneous and combination therapy



### COMBINATION THERAPY DEVICE

- 2 independent channels
- All currents
- 7" Touch Screen
- Continuous and pulsed multi-frequency treatment head, 1&3MHz, 4cm<sup>2</sup>

## Gymna 400 Series

Combination device for 2-channel electrotherapy, 2-pole and 4-pole.  
Diagnostics, ultrasound, laser, simultaneous and combination therapy



### COMBINATION THERAPY DEVICE

- 2 independent channels
- All Treatment currents
- 10.4" Full Color Touch Screen
- Anatomical library
- Mobile integrated trolley
- Continuous and pulsed multi-frequency treatment head, 1&3MHz, 4cm<sup>2</sup>



## Gymna ShockMaster 500

### RADIAL SHOCKWAVE DEVICE

- Silent compressor
- 7" Full color LCD
- Frequency: 0,5-21Hz.(35 Vector)
- Pressure:0,3-5 bar
- Anatomical Library (GTS)
- Mobile Trolley
- 2.5Millions shock/Revision



## Gymna Cryoflow ICE-CT

- **35 pre-programmed settings** based on indications or objectives
- **User friendly:** no refilling
- **Short treatment duration** due to powerful temperature reduction
- **Various nozzles** depending on the area to be treated
- **2 treatment methods:** static and handsfree or dynamic and manual



## Gymna Go

Multi sections treatment table with backrest and electric height adjustment.

### Specifications:

- Dimensions : 205 x 67 cm
- Height adjustment : 50 – 100 cm
- Lifting capacity: 200 kg



## Gymna W-Move

The Gymna W-Move is a total training concept of nine devices, providing different types of exercises that will increase your patients' flexibility, physical strength and well-being. This type of circuit training allows many people to participate at the same time and you don't even need a lot of space



## Gymna ShockMaster 300



Ergonomically shaped, lightweight handpiece

Full colour 7 inch touchscreen for intuitive operation

Optional trolley with two shelves, a drawer and integrated bottle holder



## M8

Robotized optical group with 3 MLS® sources with power up to 3.6 W – Peak Power 3 x 75 W  
Handpiece with optical group consisting of MLS® source with power up to 1.2 W – Peak Power 75 W  
Automatic calculation of the emitted energy according to the set parameters

4 treatment modalities:

- A) pain management
- B) trigger points
- C) anti-edema
- D) biostimulation



## Mphi 75

Handpiece with 1 MLS® source – 75 W Peak Power Ø2 cm,  
Target Area 3 cm<sup>2</sup>  
Automatic calculation of the emitted energy according to the set parameters


4 treatment modalities:

- A) pain management
- B) trigger points
- C) anti-edema
- D) biostimulation



## SHI

Hilterapia® is a therapeutic technique based on an exclusive high-intensity laser emission (HILT® = High Intensity Laser Therapy) that, working in absolute safety, allows to obtain an effective in depth therapeutic action.



## HIRO TT

HIGH POWER LASER + COOLING

- High Intensity pulsed Nd:YAG
- Peak Power (max): 3,000 W
- Intensity (max): 15,000 W/cm<sup>2</sup>
- Average power (max): 10.5 W
- Auto-calibration of the laser
- Smart Cooler Cooling system: 18C° to 28C°



## ASA Magnetotherapy

MAGNETO THERAPY DEVICE (FULL BODY AUTO SCAN)

- Frequency from 0.5 to 100 Hz
- Magnetic field intensity variable from 5 to 100%
- Treatment time from 1 to 99'
- Pre-set, Adjustable, Saveable programs.
- Electromechanical system for solenoid movement managed by microprocessor
- Motorized couch with Ø 80 cm solenoid diameter



## EASY Qs

PORTABLE AND CUSTOMISABLE. EASY Qs is ideal for the specific treatment of localised areas thanks to theFlexa applicator.



## SYSTEM 4 PRO™

### MULTI JOINT SYSTEM-ISOKIETEC

#### Modes:

Isokinetic Resistance, Isometric, Isotonic  
Passive Motion, Reactive Eccentric.

- Concentric speed: 500 deg/sec
- Eccentric speed: 300 deg/sec
- Concentric torque: 500 ft-lb (680 Nm)
- Eccentric torque: 400 ft-lb (542 Nm)
- Passive speed: 0.25 deg/sec
- Passive torque: 0.5 ft-lb (0.68 Nm)
- Isotonic torque: 0.5 ft-lb (0.68 Nm)
- Chair Height: Motorized adjustment

#### Attachments:

Knee, Shoulder, Ankle, Elbow, Wrist, Hip



## Balance System SD

Stability Levels: 12 dynamic levels, plus  
locked for static measurements

Platform Tilt: 20° from horizontal in all  
directions.

Support Handles:

Adjustable from

26" - 36" above platform (66 - 91 cm)

User Capacity: 60 - 400 lb (27 - 182 kg)



## GAIT TRAINING + UNWEIGHING SYSTEM.

The Biodesx NxStep™ Unweighing System enables partial weight-bearing therapy with open access to the patient. Offloading a percentage of body weight offers the opportunity for early rehabilitation while providing a safe environment for patient and therapist. The only treadmill with an instrumented deck that monitors and records step length, stepspeed and step symmetry, the Gait Trainer 3 lets you prove patients are getting better, faster – and documents results.




## Lokomat®Pro

Robot-assisted therapy enables effective and intensive training and ensures the optimal exploitation of neuroplasticity and recovery potential.



## C-Mill®VR+

C-Mill is a unique development in which a theoretical innovation has led to a practical immediately applicable and user-friendly product for rehabilitation.



## Erigo®Pro

The Erigo enables a single therapist to provide mobilization, verticalization, and sensorimotor stimulation at the same time, safely and efficiently.



## Armeo®

The Armeo devices improve therapy efficiency by reducing the therapist's physical effort and the need for continuous therapeutic guidance.





## ACTIV CT9

The CT9 device converts electrical current into a fixed stable frequency of 448 kHz, increasing ionic exchange at a cellular level which encourages tissue repair and assists in pain management.

## ACTIV 701

The use of INDIBA® Activ 701 improves recovery and reduces pain. The body's natural repair process will be accelerated and treatment of an injury can begin in the acute phase; reducing pain, swelling and speeding up recovery.

## TECHNICAL SPECIFICATIONS

	ACTIV 701	ACTIV 801	ACTIV CT8	ACTIV CT9
Output frequency	448 kHz	448 kHz	448kHz ±1kHz	448 kHz ±1kHz
Max. output power in RES mode	65 W	100 W	100 W	200 W
Max. output power in CAP mode	250 VA	350 VA	350 VA	450 VA
Continuous modulation	-	-	20 kHz	20 kHz
Return plate warmer	-	-	Integrated in the column	
Pre-set protocols	-	-	YES	
Customised protocols	-	-	YES	
Monitoring software	-	-	Optional	YES



## LYMPHATIC DRAINAGE 8 CH

### Features:

- 8 Chamber System
- Gradient, sequential pressure
- Distal to proximal movement of fluid
- Unilateral or bilateral use
- Slim and light design
- Very quiet
- Compliance meter
- Individual pressure adjustment in each chamber



## The PiezoWave2

Piezo shockwave therapy sources are available with single-layer or double-layer technology, depending on the required energy levels, with peak pressures of up to 82 MPa.



### Focused shockwave

The focused therapy source of the PiezoWave<sup>2</sup> is characterized by its well defined and precise focal zone and a penetration depth of up to 100 mm (PiezoWave<sup>2</sup> FB10 G10 therapy source). Its guaranteed life span of 5 million pulses is unsurpassed.



### Linear focused shockwave

Breaking new ground in shockwave therapy, the linear focused therapy source is a double-layer piezoelectric source with a unique geometry and new features. The focal zone is of a 3-dimensional area with an outstanding size of 46 mm x 20 mm x 4 mm. The energy output achieved with this therapy source has not been previously achieved and sets a new standard: called Volume ESWT.



**KNEE CPM™**

An unbeatable combination of great design and latest USB technology


- Intuitive remote control
- Unique features
- ROM from -10° to 120°
- KINETEC DATA CAPTURE™ USB technology ready



**SHOULDER CPM™**

The Anatomical Shoulder CPM Machine designed for use in professional environments

- Complete range of movements available.
- Synchronized and Anatomical motion
- Evolutive, build-to-your-need device with up to 3 extra modules available



**UPPER&LOWER ACTIVE/PASSIVE TRAINER™**

The intelligent exercise therapy system for upper and lower limbs, featuring a wide range of exercise modes and benefits for the user.



**Anti-Gravity Treadmill®**

The AlterG Anti-Gravity Treadmill® creates a powerful lifting force using Advanced Differential Air Pressure Technology (ADAPT) that allows for ultra-low impact ambulation. The user wears neoprene shorts and zips into a pressurized, airtight enclosure suspended over the treadmill surface.



## CARDIOVIT CS-200 OFFICE ERGOSPIRO

Covers all Cardiopulmonary Exercise Test applications. Breath by breath gas exchange analysis (VE, VO2, VCO2) with unique maintenance free ultrasonic cell for CO2 measurement Fully integrated 12-lead stress ECG with wireless acquisition improves patient comfort Full disclosure storage of unfiltered, continuous ECG signals.



## Manuthera 242

All-in-one treatment table for demanding examination and treatment use by various manual therapists, such as sports massage therapists, naprapaths, osteopaths, OMT and physical therapists and chiropractors.



## HOTPACKS UNIT

Whitehall Thermalators have the best elements of a hot therapy unit to maintain a constant water temperature for comfortable effective therapy.



## Hubbard Tanks

The 425-gallon Hubbard Tank is most commonly used to treat the entire body simultaneously— arms, legs, and trunk of the body. But unlike regular hydrotherapy whirlpools, Hubbard tanks' figure-eight shape gives medical staff easier access to patients.



USA

## Electra

The limb galvanic bath Electra is designed to provide application of galvanic current in water on the upper and lower extremities. While flowing from one limb to another, the galvanic current causes polarization of the tissues.



## HUBBARD

A wholebody hydrotherapy bath tub with manual controls. It is suitable for hydromassage and subaquatic massage treatments. Its special shape and large volume make it suitable for individual patient exercise.



SLOVAKIA

**gymna**<sup>®</sup>

**AST**  
A S A L A S E R  
Research and Therapeutic Solutions

**BIODEX**  
A MIRION MEDICAL COMPANY

**Hocoma**

**INDIBA**<sup>®</sup>  
ACTIV THERAPY  
MEDICAL PROIONIC SYSTEM

WHITEHALL  
SINCE 1946  
SPORTS MEDICINE

**CHIRANA PROGRESS**  
MEDIPROGRESS GROUP

**SCHILLER**  
The Art of Diagnostics

**ALTER** **G**  
Anti-Gravity Treadmill<sup>®</sup>

**BIO** **BIO COMPRESSION SYSTEMS INC.**  
#1 in Compression Therapy

**kinetec**<sup>®</sup>

**LOJER**<sup>®</sup>

**RICHARD WOLF**  
spirit of excellence

**ELvation**<sup>®</sup>  
think in value

**FEI**  
FABRICATION ENTERPRISES INC.

**FP**  
PERFORMANCE HEALTH<sup>™</sup>



**WWW.LMC-EG.COM**

